## **▶** Virus purification and detection products

## **Pannarase Nucleases**

Effective viscosity reduction and removal of nucleic acids from protein solutions

### **Product Introduction**

Pannarase Nuclease is a genetically engineered endonuclease from *Serratia marcescens*, which can completely digests nucleic acids to 5′—monophosphate teminated oligonucleotides 2 to 5 bases in length. It degrades all forms of DNA and RNA (single stranded, double stranded, linear and circular) while having no proteolytic activity.

The enzyme is functional between pH 6 and 10 and from 0–42°C. The ability of Pannarase to rapidly hydrolyze nucleic acids makes the enzyme an excellent choice for viscosity reduction to reduce processing time and increase yields of protein.

#### **Product Advantages**

- Animal origin—free, without ampicillin
- High purity(>99%),Excellent unit enzyme activity
- Advanced manufacturing process, non-traditional His-tag purification, exclusion of metal ion risks
- Low endotoxin level and strict quality control criteria to ensure unit enzyme activity and batch—to—batch stability


### **Product Application**

Widely used in virus purification processes, such as lentivirus, adenovirus, lysovirus, etc. Removal of nucleic acids from biological products. Reduce the viscosity of cell lysates and cell supernatants to saving processing time and improve protein yields


# SEC-HPLC purity

determination


**Test Data** 



Pannarase SEC-HPLC purity analysis SEC-HPLC reault showed no aggregate formation, and the purity reached 99.72%.



As analyzed by LC-MS, the protein molecule was consistent with the theoretical value, there was no impurity peak, and the structure was homogeneous.



The plasmid DNA concentration was 1 mg/ml in the 20 uL reaction system. Compared with competitors, Pannarase Nuclease had excellent digestion activity.

#### **Order Options**

Cat No.	Grade	Name	Size
CYG002F0010	GMP Grade	Pannarase Nucleases	100,000 U
CYG002F0050	GMP Grade	Pannarase Nucleases	500,000 U
CYG002F0500	GMP Grade	Pannarase Nucleases	5,000,000 U
CY002F0010	RUO(Research Use Only) Grade	Pannarase Nucleases	100,000 U
CY002F0050	RUO(Research Use Only) Grade	Pannarase Nucleases	500,000 U
CY002F0500	RUO(Research Use Only) Grade	Pannarase Nucleases	5,000,000 U